



MAP3K13 Antibody, Biotin conjugated

Product Code	CSB-PA013421LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43283
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 13 protein (506-835AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 13 (EC 2.7.11.25) (Leucine zipper-bearing kinase) (Mixed lineage kinase) (MLK), MAP3K13, LZK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP3K13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.