



MAPK6 Antibody, Biotin conjugated

Product Code	CSB-PA013458LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16659
Immunogen	Recombinant Human Mitogen-activated protein kinase 6 protein (383-712AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase 6 (MAP kinase 6) (MAPK 6) (EC 2.7.11.24) (Extracellular signal-regulated kinase 3) (ERK-3) (MAP kinase isoform p97) (p97-MAPK), MAPK6, ERK3 PRKM6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAPK6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.