



MARK3 Antibody, HRP conjugated

Product Code	CSB-PA013498LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27448
Immunogen	Recombinant Human MAP/microtubule affinity-regulating kinase 3 protein (368-667AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MAP/microtubule affinity-regulating kinase 3 (EC 2.7.11.1) (C-TAK1) (cTAK1) (Cdc25C-associated protein kinase 1) (ELKL motif kinase 2) (EMK-2) (Protein kinase STK10) (Ser/Thr protein kinase PAR-1) (Par-1a) (Serine/threonine-protein kinase p78), MARK3, CTAK1 EMK2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MARK3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.