



MATN4 Antibody

Product Code	CSB-PA013523LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95460
Immunogen	Recombinant Human Matrilin-4 protein (326-506AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Matrilin-4, MATN4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MATN4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.