



MBP Antibody, HRP conjugated

Product Code	CSB-PA013551LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02686
Immunogen	Recombinant Human Myelin basic protein (1-197AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myelin basic protein (MBP) (Myelin A1 protein) (Myelin membrane encephalitogenic protein), MBP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MBP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.