



MGST3 Antibody, Biotin conjugated

Product Code	CSB-PA013793LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14880
Immunogen	Recombinant Human Microsomal glutathione S-transferase 3 protein (30-70AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Microsomal glutathione S-transferase 3 (Microsomal GST-3) (EC 2.5.1.18) (Microsomal GST-III), MGST3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MGST3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.