



MNDA Antibody, HRP conjugated

Product Code	CSB-PA014688LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41218
Immunogen	Recombinant Human Myeloid cell nuclear differentiation antigen protein (84-233AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myeloid cell nuclear differentiation antigen, MNDA
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	MNDA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.