



MYF6 Antibody, Biotin conjugated

Product Code	CSB-PA015288LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23409
Immunogen	Recombinant Human Myogenic factor 6 protein (26-232AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myogenic factor 6 (Myf-6) (Class C basic helix-loop-helix protein 4) (bHLHc4) (Muscle-specific regulatory factor 4), MYF6, BHLHC4 MRF4
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MYF6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.