



NAMPT Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA015422EB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P43490 |
| Immunogen | Recombinant Human Nicotinamide phosphoribosyltransferase protein (1-280AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Nicotinamide phosphoribosyltransferase (NAMPRase) (Nampt) (EC 2.4.2.12) (Pre-B-cell colony-enhancing factor 1) (Pre-B cell-enhancing factor) (Visfatin), NAMPT, PBEF PBEF1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cardiovascular |
| Target Names | NAMPT |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |