



NAMPT Antibody, Biotin conjugated

Product Code	CSB-PA015422ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43490
Immunogen	Recombinant Human Nicotinamide phosphoribosyltransferase protein (1-280AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nicotinamide phosphoribosyltransferase (NAMPRase) (Nampt) (EC 2.4.2.12) (Pre-B-cell colony-enhancing factor 1) (Pre-B cell-enhancing factor) (Visfatin), NAMPT, PBEF PBEF1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	NAMPT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.