



NAT1 Antibody, Biotin conjugated

Product Code	CSB-PA015464LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18440
Immunogen	Recombinant Human Arylamine N-acetyltransferase 1 protein (1-288AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Arylamine N-acetyltransferase 1 (EC 2.3.1.5) (Arylamide acetylase 1) (Monomorphic arylamine N-acetyltransferase) (MNAT) (N-acetyltransferase type 1) (NAT-1), NAT1, AAC1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NAT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.