



NOS1AP Antibody, Biotin conjugated

Product Code	CSB-PA015942LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75052
Immunogen	Recombinant Human Carboxyl-terminal PDZ ligand of neuronal nitric oxide synthase protein (176-288AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carboxyl-terminal PDZ ligand of neuronal nitric oxide synthase protein (C-terminal PDZ ligand of neuronal nitric oxide synthase protein) (Nitric oxide synthase 1 adaptor protein), NOS1AP, CAPON KIAA0464
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NOS1AP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.