



NR0B1 Antibody, Biotin conjugated

Product Code	CSB-PA016041LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51843
Immunogen	Recombinant Human Nuclear receptor subfamily 0 group B member 1 protein (1-243AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear receptor subfamily 0 group B member 1 (DSS-AHC critical region on the X chromosome protein 1) (Nuclear receptor DAX-1), NR0B1, AHC DAX1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NR0B1