



NUCB2 Antibody, Biotin conjugated

Product Code	CSB-PA016146HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80303
Immunogen	Recombinant Human Nucleobindin-2 protein (25-106AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleobindin-2 (DNA-binding protein NEFA) (Gastric cancer antigen Zg4) (Prepronesfatin) [Cleaved into: Nesfatin-1], NUCB2, NEFA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NUCB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.