



PCCA Antibody, Biotin conjugated

Product Code	CSB-PA017522LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05165
Immunogen	Recombinant Human Propionyl-CoA carboxylase alpha chain, mitochondrial protein (585-720AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Propionyl-CoA carboxylase alpha chain, mitochondrial (PCCase subunit alpha) (EC 6.4.1.3) (Propanoyl-CoA:carbon dioxide ligase subunit alpha), PCCA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PCCA