



PCCA Antibody, Biotin conjugated

Product Code	CSB-PA017522LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05165
Immunogen	Recombinant Human Propionyl-CoA carboxylase alpha chain, mitochondrial protein (585-720AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Propionyl-CoA carboxylase alpha chain, mitochondrial (PCCase subunit alpha) (EC 6.4.1.3) (Propanoyl-CoA:carbon dioxide ligase subunit alpha), PCCA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PCCA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.