



PCGF2 Antibody, HRP conjugated

Product Code	CSB-PA017606LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35227
Immunogen	Recombinant Human Polycomb group RING finger protein 2 protein (103-316AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polycomb group RING finger protein 2 (DNA-binding protein Mel-18) (RING finger protein 110) (Zinc finger protein 144), PCGF2, MEL18 RNF110 ZNF144
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PCGF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.