



Pcsk1 Antibody

Product Code	CSB-PA017640LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28840
Immunogen	Recombinant Rat Neuroendocrine convertase 1 protein (111-752AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuroendocrine convertase 1 (NEC 1) (EC 3.4.21.93) (Prohormone convertase 1) (Proprotein convertase 1) (PC1), Pcsk1, Bdp Nec-1 Nec1
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Pcsk1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.