



PHF21A Antibody, HRP conjugated

Product Code	CSB-PA017911LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96BD5
Immunogen	Recombinant Human PHD finger protein 21A protein (312-487AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PHD finger protein 21A (BHC80a) (BRAF35-HDAC complex protein BHC80), PHF21A, BHC80 KIAA1696
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PHF21A