



PHF21A Antibody, Biotin conjugated

Product Code	CSB-PA017911LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96BD5
Immunogen	Recombinant Human PHD finger protein 21A protein (312-487AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PHD finger protein 21A (BHC80a) (BRAAF35-HDAC complex protein BHC80), PHF21A, BHC80 KIAA1696
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PHF21A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.