



PIK3R2 Antibody, FITC conjugated

Product Code	CSB-PA018005LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00459
Immunogen	Recombinant Human Phosphatidylinositol 3-kinase regulatory subunit beta protein (75-125AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylinositol 3-kinase regulatory subunit beta (PI3-kinase regulatory subunit beta) (PI3K regulatory subunit beta) (PtdIns-3-kinase regulatory subunit beta) (Phosphatidylinositol 3-kinase 85 kDa regulatory subunit beta) (PI3-kinase subunit p85-beta) (PtdIns-3-kinase regulatory subunit p85-beta), PIK3R2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PIK3R2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.