



# PIK3R2 Antibody, Biotin conjugated

|                     |                                                                                                                                                                                                                                                                                                                        |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code        | CSB-PA018005LD01HU                                                                                                                                                                                                                                                                                                     |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                                                                                                                                          |
| Uniprot No.         | O00459                                                                                                                                                                                                                                                                                                                 |
| Immunogen           | Recombinant Human Phosphatidylinositol 3-kinase regulatory subunit beta protein (75-125AA)                                                                                                                                                                                                                             |
| Raised In           | Rabbit                                                                                                                                                                                                                                                                                                                 |
| Species Reactivity  | Human                                                                                                                                                                                                                                                                                                                  |
| Tested Applications | ELISA                                                                                                                                                                                                                                                                                                                  |
| Form                | Liquid                                                                                                                                                                                                                                                                                                                 |
| Conjugate           | Biotin                                                                                                                                                                                                                                                                                                                 |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                                                                                                                                                                                                                                       |
| Purification Method | >95%, Protein G purified                                                                                                                                                                                                                                                                                               |
| Isotype             | IgG                                                                                                                                                                                                                                                                                                                    |
| Clonality           | Polyclonal                                                                                                                                                                                                                                                                                                             |
| Alias               | Phosphatidylinositol 3-kinase regulatory subunit beta (PI3-kinase regulatory subunit beta) (PI3K regulatory subunit beta) (PtdIns-3-kinase regulatory subunit beta) (Phosphatidylinositol 3-kinase 85 kDa regulatory subunit beta) (PI3-kinase subunit p85-beta) (PtdIns-3-kinase regulatory subunit p85-beta), PIK3R2 |
| Immunogen Species   | Homo sapiens (Human)                                                                                                                                                                                                                                                                                                   |
| Research Area       | Signal Transduction                                                                                                                                                                                                                                                                                                    |
| Target Names        | PIK3R2                                                                                                                                                                                                                                                                                                                 |
| Usage               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                                                                                                                                                                                                                                            |