



PLEK Antibody

Product Code	CSB-PA018152LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08567
Immunogen	Recombinant Human Pleckstrin protein (38-337AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pleckstrin (Platelet 47 kDa protein) (p47), PLEK, P47
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PLEK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.