



PNLIPRP2 Antibody, Biotin conjugated

Product Code	CSB-PA018264LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54317
Immunogen	Recombinant Human Pancreatic lipase-related protein 2 protein (305-460AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pancreatic lipase-related protein 2 (PL-RP2) (EC 3.1.1.26) (EC 3.1.1.3) (Galactolipase), PNLIPRP2, PLRP2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PNLIPRP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.