



PPP1CA Antibody, Biotin conjugated

Product Code	CSB-PA018503LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62136
Immunogen	Recombinant Human Serine/threonine-protein phosphatase PP1-alpha catalytic subunit protein (192-330AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein phosphatase PP1-alpha catalytic subunit (PP-1A) (EC 3.1.3.16), PPP1CA, PPP1A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPP1CA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.