



Ppp1cb Antibody, Biotin conjugated

Product Code	CSB-PA018504LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62141
Immunogen	Recombinant Mouse Serine/threonine-protein phosphatase PP1-beta catalytic subunit protein (2-327AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein phosphatase PP1-beta catalytic subunit (PP-1B) (EC 3.1.3.16) (EC 3.1.3.53), Ppp1cb
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Ppp1cb
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.