



PPP6C Antibody, Biotin conjugated

Product Code	CSB-PA018583LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00743
Immunogen	Recombinant Human Serine/threonine-protein phosphatase 6 catalytic subunit protein (40-211AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein phosphatase 6 catalytic subunit (PP6C) (EC 3.1.3.16) [Cleaved into: Serine/threonine-protein phosphatase 6 catalytic subunit, N-terminally processed], PPP6C, PPP6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPP6C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.