



PPT1 Antibody, FITC conjugated

Product Code	CSB-PA018587LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50897
Immunogen	Recombinant Human Palmitoyl-protein thioesterase 1 protein (100-258AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Palmitoyl-protein thioesterase 1 (PPT-1) (EC 3.1.2.22) (Palmitoyl-protein hydrolase 1), PPT1, CLN1 PPT
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.