



PPT1 Antibody, Biotin conjugated

Product Code	CSB-PA018587LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50897
Immunogen	Recombinant Human Palmitoyl-protein thioesterase 1 protein (100-258AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Palmitoyl-protein thioesterase 1 (PPT-1) (EC 3.1.2.22) (Palmitoyl-protein hydrolase 1), PPT1, CLN1 PPT
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPT1