



PRKACB Antibody, FITC conjugated

Product Code	CSB-PA018689LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22694
Immunogen	Recombinant Human cAMP-dependent protein kinase catalytic subunit beta protein (2-53AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	cAMP-dependent protein kinase catalytic subunit beta (PKA C-beta) (EC 2.7.11.11), PRKACB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PRKACB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.