



PRKAR2A Antibody, HRP conjugated

Product Code	CSB-PA018697LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13861
Immunogen	Recombinant Human cAMP-dependent protein kinase type II-alpha regulatory subunit protein (43-205AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	cAMP-dependent protein kinase type II-alpha regulatory subunit, PRKAR2A, PKR2 PRKAR2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PRKAR2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.