



PRKAR2A Antibody, Biotin conjugated

Product Code	CSB-PA018697LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13861
Immunogen	Recombinant Human cAMP-dependent protein kinase type II-alpha regulatory subunit protein (43-205AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	cAMP-dependent protein kinase type II-alpha regulatory subunit, PRKAR2A, PKR2 PRKAR2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PRKAR2A