



PRKD3 Antibody, FITC conjugated

Product Code	CSB-PA018713LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94806
Immunogen	Recombinant Human Serine/threonine-protein kinase D3 protein (326-510AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase D3 (EC 2.7.11.13) (Protein kinase C nu type) (Protein kinase EPK2) (nPKC-nu), PRKD3, EPK2 PRKCN
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PRKD3