




Tmprss15 Antibody

Product Code	CSB-PA018824LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P97435
Immunogen	Recombinant Mouse Enteropeptidase protein (830-1069AA)
Raised In	Rabbit
Species Reactivity	Mouse, Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Enteropeptidase (EC 3.4.21.9) (Enterokinase) (Serine protease 7) (Transmembrane protease serine 15) [Cleaved into: Enteropeptidase non-catalytic heavy chain; Enteropeptidase catalytic light chain], Tmprss15, Entk Prss7
Immunogen Species	Mus musculus (Mouse)
Research Area	Cell Biology
Target Names	Tmprss15

Image



Western Blot

Positive WB detected in: HepG2 whole cell lysate, Rat liver tissue

All lanes: Tmprss15 antibody at 3μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 119 kDa

Observed band size: 119 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.