




PSEN2 Antibody

Product Code	CSB-PA018847LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49810
Immunogen	Recombinant Human Presenilin-2 protein (7-77AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Presenilin-2 (PS-2) (EC 3.4.23.-) (AD3LP) (AD5) (E5-1) (STM-2) [Cleaved into: Presenilin-2 NTF subunit; Presenilin-2 CTF subunit], PSEN2, AD4 PS2 PSNL2 STM2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PSEN2

Image



Western Blot

Positive WB detected in: Rat spleen tissue, Rat lung tissue

All lanes: PSEN2 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 51, 47 kDa

Observed band size: 51 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.