



PSMA1 Antibody, FITC conjugated

Product Code	CSB-PA018865LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25786
Immunogen	Recombinant Human Proteasome subunit alpha type-1 protein (163-263AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proteasome subunit alpha type-1 (EC 3.4.25.1) (30 kDa prosomal protein) (PROS-30) (Macropain subunit C2) (Multicatalytic endopeptidase complex subunit C2) (Proteasome component C2) (Proteasome nu chain), PSMA1, HC2 NU PROS30 PSC2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.