



PSMA5 Antibody, FITC conjugated

Product Code	CSB-PA018870LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28066
Immunogen	Recombinant Human Proteasome subunit alpha type-5 protein (124-180AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proteasome subunit alpha type-5 (EC 3.4.25.1) (Macropain zeta chain) (Multicatalytic endopeptidase complex zeta chain) (Proteasome zeta chain), PSMA5
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMA5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.