



PSMC4 Antibody, FITC conjugated

Product Code	CSB-PA018894LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43686
Immunogen	Recombinant Human 26S proteasome regulatory subunit 6B protein (112-181AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	26S proteasome regulatory subunit 6B (26S proteasome AAA-ATPase subunit RPT3) (MB67-interacting protein) (MIP224) (Proteasome 26S subunit ATPase 4) (Tat-binding protein 7) (TBP-7), PSMC4, MIP224 TBP7
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMC4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.