



PSMC5 Antibody, HRP conjugated

Product Code	CSB-PA018895LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62195
Immunogen	Recombinant Human 26S proteasome regulatory subunit 8 protein (253-324AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	26S proteasome regulatory subunit 8 (26S proteasome AAA-ATPase subunit RPT6) (Proteasome 26S subunit ATPase 5) (Proteasome subunit p45) (Thyroid hormone receptor-interacting protein 1) (TRIP1) (p45/SUG), PSMC5, SUG1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMC5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.