



# PTF1A Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA018967LB01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q7RTS3
<b>Immunogen</b>	Recombinant Human Pancreas transcription factor 1 subunit alpha protein (224-325AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Pancreas transcription factor 1 subunit alpha (Class A basic helix-loop-helix protein 29) (bHLHa29) (Pancreas-specific transcription factor 1a) (bHLH transcription factor p48) (p48 DNA-binding subunit of transcription factor PTF1) (PTF1-p48), PTF1A, BHLHA29 PTF1P48
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Neuroscience
<b>Target Names</b>	PTF1A
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.