



PTF1A Antibody, Biotin conjugated

Product Code	CSB-PA018967LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7RTS3
Immunogen	Recombinant Human Pancreas transcription factor 1 subunit alpha protein (224-325AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pancreas transcription factor 1 subunit alpha (Class A basic helix-loop-helix protein 29) (bHLHa29) (Pancreas-specific transcription factor 1a) (bHLH transcription factor p48) (p48 DNA-binding subunit of transcription factor PTF1) (PTF1-p48), PTF1A, BHLHA29 PTF1P48
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PTF1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.