



RAB27A Antibody, FITC conjugated

Product Code	CSB-PA019177LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51159
Immunogen	Recombinant Human Ras-related protein Rab-27A protein (2-221AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein Rab-27A (Rab-27) (GTP-binding protein Ram), RAB27A, RAB27
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RAB27A