



ROCK2 Antibody, HRP conjugated

Product Code	CSB-PA020059DB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75116
Immunogen	Recombinant Human Rho-associated protein kinase 2 protein (1109-1388AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho-associated protein kinase 2 (EC 2.7.11.1) (Rho kinase 2) (Rho-associated, coiled-coil-containing protein kinase 2) (Rho-associated, coiled-coil-containing protein kinase II) (ROCK-II) (p164 ROCK-2), ROCK2, KIAA0619
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ROCK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.