



RPL23A Antibody, HRP conjugated

Product Code	CSB-PA020187LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62750
Immunogen	Recombinant Human 60S ribosomal protein L23a protein (1-92AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	60S ribosomal protein L23a (Large ribosomal subunit protein uL23), RPL23A
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RPL23A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.