



RPS4Y1 Antibody, HRP conjugated

Product Code	CSB-PA020457LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22090
Immunogen	Recombinant Human 40S ribosomal protein S4, Y isoform 1 protein (197-263AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	40S ribosomal protein S4, Y isoform 1 (Small ribosomal subunit protein eS4), RPS4Y1, RPS4Y
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RPS4Y1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.