





SERPINC1 Antibody

Product Code	CSB-PA021079LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01008
Immunogen	Recombinant Human Antithrombin-III protein (37-180AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:3000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Antithrombin-III (ATIII) (Serpine C1), SERPINC1, AT3
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SERPINC1


Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA021079LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA021079LA01HU at dilution of 1:100

**Western Blot**

Positive WB detected in: A549 whole cell lysate, Raji whole cell lysate, MCF-7 whole cell lysate, Rat spleen tissue

All lanes: SERPINC1 antibody at 2.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 53 kDa

Observed band size: 53 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.