



SERPINC1 Antibody, Biotin conjugated

Product Code	CSB-PA021079LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01008
Immunogen	Recombinant Human Antithrombin-III protein (37-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Antithrombin-III (ATIII) (Serp C1), SERPINC1, AT3
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SERPINC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.