



SFRP2 Antibody, FITC conjugated

Product Code	CSB-PA021139LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96HF1
Immunogen	Recombinant Human Secreted frizzled-related protein 2 protein (68-186AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Secreted frizzled-related protein 2 (FRP-2) (sFRP-2) (Secreted apoptosis-related protein 1) (SARP-1), SFRP2, FRP2 SARP1
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	SFRP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.