



SLC16A8 Antibody, Biotin conjugated

Product Code	CSB-PA021415LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95907
Immunogen	Recombinant Human Monocarboxylate transporter 3 protein (432-504AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Monocarboxylate transporter 3 (MCT 3) (Solute carrier family 16 member 8), SLC16A8, MCT3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC16A8