



SLC5A1 Antibody, FITC conjugated

Product Code	CSB-PA021675LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13866
Immunogen	Recombinant Human Sodium/glucose cotransporter 1 protein (564-643AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium/glucose cotransporter 1 (Na(+)/glucose cotransporter 1) (High affinity sodium-glucose cotransporter) (Solute carrier family 5 member 1), SLC5A1, NAGT SGLT1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC5A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.